

# PROTEC BATHS

CRAFTED IN BRITAIN

## Over Bath Stretcher Technical Information


**The Protec Baths Ltd Over Bath Stretcher is designed to be used as a personal support platform to be placed on the top of an assisted bath manufactured by Protec Baths Ltd.**

It is easy and quick to deploy, making assisting the user easier for carers.

With the ability to be used on an assisted bath with and without a power traverse seat, it means the bath can be made accessible for a wider range of different users.


### **Benefits of using the Over Bath Stretcher:**

- Makes the Bathing experience more accessible for more types of users
- Saves the need for additional changing tables thus saving money
- Ability to dry or tend to the user within the bath perimeter saving much needed space in the bathroom
- Potential to use the stretcher as a shower trolley to wash user over the bath



# PROTEC BATHS

CRAFTED IN BRITAIN



Overall Size – 1700 x 700mm

A – 600mm

B – 1520mm

## Features:

- **Removable cover for washing or replacement over time**
- **Antimicrobial & FR Fabric**
- **Drainage holes in length of the stretcher**
- **Strong tubular steel frame with supports for strength**

# PROTEC BATHS

CRAFTED IN BRITAIN

## **Cleaning Information:**

**The Frame and cover can be cleaned using mild cleaning agents and anti-bacterial solutions.**

**Any abrasive cleaners, or chemical agents and solvents must not be used, for example, Drain/oven Cleaner, Acetone, or Acids.**

**Please always pay attention to instruction labels on the outside of the cleaning agent containers.**